

Minncare® Test Strips

Validated Peracetic Acid Test Strips

Minncare Test Strips provide quick results with easy-to-read indicators. Use the Minncare 1% TS Test Strips as a pass/fail measurement for adequate concentration of sterilant after dilution. Then verify residual levels after rinse-out using Minncare Residual Test Strips.

Minncare Peracetic Acid Test Strips represent the leading edge of test strip technology. Designed by the producers of Minncare Cold Sterilant, these test strips combine ease-of-use with precise results.

Both the Minncare 1% TS and the Minncare Residual Test Strips may be used up to and including the last day of the month indicated by the expiration printed on the container even after opening. There is no need to separately record and verify container opening dates.

Minncare® 1% Test Strips

Minncare 1% Test Strips verify adequate concentration of a 1% Minncare solution during sanitization of your high purity water system. These convenient dip-and-read strips can be used at any test port in your system to provide an accurate pass/fail measurement of Minncare solution.



Minncare® Residual Test Strips

Minncare Residual Test Strips are easy to use and are conveniently packaged 100 test strips per container. Each strip has a “reaction zone” pad located at one end of the strip. When stored and used properly, the test strip pad can indicate residual Minncare levels as low as 1 part per million (ppm). The reaction zones are set to provide six different reading levels at 100, 30, 10, 3, 1 and 0 (zero) parts per million.

| Minncare Ordering Information | | | | | |
|-------------------------------|-------------------|-------------------|-----------------------|---------------------------------------|----------------------------|
| U.S. Part No. | Canadian Part No. | European Part No. | Asia Pacific Part No. | Description | Packaging |
| 176-01-001 | 78399-647 | n/a | n/a | Minncare Cold Sterilant 4 x 1 quart | 4 x 1 quarts |
| 176-01-002 | 78399-646 | n/a | n/a | Minncare Cold Sterilant 4 US gallons | 4 x 1 US gallons |
| 176-01-003 | 78399-648 | n/a | n/a | Minncare Cold Sterilant 40 US gallons | 40 gallons in 55 gal. drum |
| n/a | n/a | 176-01-125 | 78397-983 | Minncare Cold Sterilant 2 x 5 liters | 2 x 5 liters |
| n/a | n/a | 176-01-121 | 78397-825 | Minncare Cold Sterilant 2 x 10 liters | 2 x 10 liters |
| n/a | n/a | 176-01-129 | 78397-229 | Minncare Cold Sterilant 6 x 1 kg | 6 x 1 kg |
| 185-40-004 | 78338-000 | 185-40-004 | n/a | Minncare Residual Test Strips | 100 per vial |
| 185-40-005 | 78339-000 | 185-40-005 | n/a | Minncare 1% Test Strips | 100 per vial |

Minncare® is a registered trademark of Minntech Corporation, a Cantel Medical Company



Mar Cor Purification
4450 Township Line Road
Skipack, PA 19474-1429
Tel: (484) 991-0220
Toll Free: (800) 346-0365
Fax: (484) 991-0230

Mar Cor Purification
14550 28th Avenue North
Plymouth, MN 55447
Tel: (484) 991-0220
Toll Free: (800) 633-3080
Fax: (763) 210-3868

SepSol Process Solutions
A UnuZucturerfbs FepresentUtij e
7\icU[o =? UIUa Unoo =8etroit
fB* -t' ((!+ * \$ k k k "SepSol"coa
inZo4 sepsol"coa

